

# EXO-NET® PAN-EXOSOME CAPTURE

## PRODUCT INFORMATION




EXO-NET Pan-Exosome Capture is a magnetic bead-based immunoaffinity EV capture device, where a 3-dimensional antibody matrix is constructed on a carbon nanoparticle. The complement of antibodies attached to the bead has been designed to capture a wide range of extracellular vesicles (EVs, including exosomes) from different cell types (pan-EV capture). EXO-NET is designed for On-Bead Analysis of captured nanoparticles and ON-Bead lysis for downstream analysis of EV-associated molecules including proteins, oligonucleotides (mRNAs and microRNAs) and lipids. EXO-NET also can be tuned to preferentially isolate EVs from specific cell types. EXO-NET is a scalable EV isolation solution for high throughput screening.

### EXO-NET® Pan-Exosome Capture Superior EV Yield

EXO-NET Pan-Exosome Capture isolates more nanoparticles from normal human plasma\*, than other commercially available bead-based kits tested (Kit A and B below), as measured by nanoparticle tracking analysis using ZetaView). The concentration of nanoparticles remaining in the input sample after EV isolation was reduced to < 20% using EXO-NET compared to >30% with the other kits.

### Depletion of human plasma nanoparticles




\*plasma was centrifuged at 10,000 g before EV isolation: Commercially available kits (Exosome Isolation Kit Pan, human - Milteni and Plasma/Serum Exosome Purification and RNA Isolation Mini Kit, Norgen. Data are presented as the mean ± SEM (n= 3).

### EXO-NET Pan-Exosome Capture Characteristics


Speed	15 minutes from sample to EV isolation
Compatibility	With downstream EV-associated analyte analyses including RT-qPCR, Western Blot, ELISA, LC/MS and NGS.
Purity and Specificity	Proprietary 3D affinity matrix for improved purity and specificity.
Scalability	Customizability and scalability for isolation of enriched exosomes
Versatile	For use with any biofluid sample

### EXO-NET® Pan-exosome – Simply Workflow




### High Purity EV Protein

EXO-NET captured EVs from pooled normal human plasma showed present of EV known protein markers CD63, CD81, CD9, Flotillin-1, and TSG101.




### Reduced Apolipoprotein contamination

EXO-NET not only outperformed other EV isolation kits regarding the enrichment of EV associated marker such as Flotillin-1 (Panel A). EXO-NET EV lysates have less co-isolated contaminants (e.g., ApoB, Panel B) when compared to other EV isolation kits (lanes #1, #2 and #3). Source: INOVIQ Collaborator Data





### Reduced Contamination by Plasma Proteins

Mass spectrometry analysis of plasma EV-lysates indicates that EXO-NET reduces albumin contamination by more than half when compared to the other commercially available kits tested. Source: INOVIQ Collaborator Data



### High Yield EV miRNA and mRNA

Overview of miRNAs (upper panel A) and mRNAs (lower panel B) yield and recovery using EXO-NET Pan-exosome and 4 other commercially available EV isolation kits (A, B, C, D) were assayed by RT-qPCR for miRNAs miR-16, let-7a and miR-21 and mRNA for GAPDH, OAZ1, RPLP0 and SERF2. EXO-NET results in equivalent or higher recovery of plasma EV RNA compared to other 4 EV isolation kits as indicated by a lower CT value.



### Product ordering information

Cat No.	Product Name	Isolations	Size
40031	EXO-NET Pan-Exosome Capture	60	1.0mL
40033	EXO-NET Pan-Exosome Capture	30	0.5mL
40036	EXO-NET Pan-Exosome Capture	15	0.25mL

Number of isolations based on 15 ul; EXO-NET per sample

### Contact Information

To learn more about EXO-NET® visit [EXO-NET](http://inoviq.com/site/products/exo-net-pan-exosome) (inoviq.com/site/products/exo-net-pan-exosome) or email [info@inoviq.com](mailto:info@inoviq.com).